




Drug dispensary bin



110% total opening shelf



Lock with key



Technical Characteristics.

- _ Aluminium and steel structure with epoxy coating
- _ ABS upper top with extended up rims on three sides
- _ Lateral push handle
- _ High impact ABS base
- _ Drawers 110% total opening, colours: red, blue, purple, white, and light green
- _ Telescopic slides with automatic shutting system
- _ Identification card holder, colour white. Available in colours red, purple, yellow, green and blue.
- _ 4 shelves with 16 drug dispensary bins 230x170x70mm, 110% total opening
- _ Drug bins with transparent dividers and identification card holder
- _ Lock with key
- _ Ø 100mm rolling wheels, two with brakes



78



46



97



60



0,34