

Travel

Big Travel

Crutch for the journey

Travel

Detachable crutch for the journey

Equipment:

- detachable – longest piece: 44 cm
- delivery inclusive bag
- 3 adjustments in height from handle to cuff
- ergonomic handle
- reflector front/backwards
- 8 adjustments in height
- compression fitting
- solid ferrule with a high diameter
- approved up to 130 kg body weight

ref.no. 108.50

BIG Travel

Double telescopable crutch for the journey

Equipment:

- 3 adjustments in height from handle to cuff
- ergonomic handle
- reflector front
- 9 adjustments in height
- double telescopable tube set
- solid ferrule with a high diameter
- approved up to 150 kg body weight

ref.no. 108.00



ref.no. 108.50

ref.no. 108.00

The travel support!

Type:	Travel	BIG Travel
height from handle to cuff:	22 cm - 25 cm	22 cm - 25 cm
height of handle:	81 cm - 95 cm	80 cm - 100 cm
diameter Ø in mm:	22/19/22	25/22/19
ref.no.:	108.50	108.00