



## duplex®

The duplex® stethoscopes from **Riester** combine first-class quality at an attractive price. They offer outstanding acoustic performance for the daily challenges in office practices and hospitals.

The acoustic performance has been considerably improved.

- Double-headed chest-piece available in two types: chromium-plated brass or specially lightweight aluminium, galvanized.
- Special membrane for perfect auscultation, Ø 48 mm (1.89 in), with particularly flat and soft non-chill rims for improved adaptation on the skin.
- Bell, Ø 36 mm (1.42 in)
- Overall length: 77 cm (30.32 in)
- duplex® is supplied with a pair of replacement ear-tips and a replacement membrane.

### Order information duplex®:

#### Aluminium chest-piece

- |             |              |
|-------------|--------------|
| No. 4001-01 | ■ black      |
| No. 4001-02 | ■ slate grey |
| No. 4031-03 | ■ blue       |
| No. 4033-05 | ■ green      |

#### Chest-piece made of chromium-plated brass

- |             |              |
|-------------|--------------|
| No. 4011-01 | ■ black      |
| No. 4011-02 | ■ slate grey |
| No. 4011-03 | ■ blue       |
| No. 4011-05 | ■ green      |

duplex®

PRODUCT CHARACTERISTICS



**duplex® baby** is a high-grade quality product from the **Riester** company. The stethoscope with double chest-piece was especially developed for babies and children. Considerably improved acoustics.

- Particularly light precision chest-piece made of aluminium.
- Special membrane, 36 mm (1.42 in), with particularly flat and soft non-chill rims for improved adaptation on the skin.
- Bell, Ø 28 mm (1.10 in)
- Overall length: 76,5 cm (30.12 in)
- **duplex® baby** is supplied with a pair of replacement ear-tips and a replacement membrane.

**Order information duplex® baby:**

No. 4041 ■ blue



**duplex® neonatal** is a quality stethoscope with precision finished double chest-piece for new-born children and infants. With considerably improved acoustics.

- Particularly small precision chest-piece made of aluminium especially for new-born children and infants.
- Special membrane, Ø 28 mm (1.10 in), with particularly flat and soft non-chill rims for improved adaptation on the skin.
- Bell, Ø 22 mm (0.87 in)
- Overall length: 76 cm (30.12 in)
- **duplex® neonatal** is supplied with a pair of replacement ear-tips and a replacement membrane.

**Order information duplex® neonatal:**

No. 4051 ■ blue

The **duplex® Teaching Stethoscope** with double chest-piece and two binaural has been developed especially for teaching purposes. It is excellently suited for auscultation by two individuals. Considerably improved acoustic performance.

- Double-headed chest-piece made of stainless steel.
- Matt chromium plated binaural with integrated spring.
- Special membrane, Ø 48 mm (1.89 in), with particularly flat and soft non-chill rims for improved adaptation on the skin.
- Bell, Ø 36 mm (1.42 in)
- Overall length: 115 cm (45.28 in)
- **duplex® Teaching Stethoscope** is supplied with two pairs of replacement ear-tips and a replacement membrane.

**Order information duplex® Teaching Stethoscope:**

No. 4002-02 ■ slate grey

